

**Grain
Flour
Food
Feed**

Quality Control & Analysis Solutions for Millers, Breeders and Growers

**19th Annual IAOM Mideast & Africa
District Conference & Expo
Arusha - Tanzania
November 17, 2008**

**Gunnar Nilsson
Area Sales Manager**

Complete Flour & Grain Analysis

- NIR – Powder & Whole Grain



- 6 Second Analysis - Diode Array NIR



- Falling Number – Sprout Damage



- RVA – Enzyme analysis, Starch Analysis



- doughLAB – Farinograph Replacement & More



- Glutomatic Gluten Tester



- Grain Moisture Testers



**Grain
Flour
Food
Feed**

Inframatic Flour Analyzer

Designed and Built for Powder and Flour Analysis:

- **Easy-to-Use routine analyzer – No PC**
- **Dust free enclosure with no cooling fan**
- **30 second accurate analysis**
 - Moisture, Protein, Ash
- **Patented Ash analysis**
- **Over 9,000 Inframatics in use today!**



IM9140 NIR Flour Analyzer

New generation NIR for mills

- **Fast, Accurate Results**
 - Moisture, Protein, Ash, Water Absorption...
- **Automated Sample Packer**
- **Remote Connectivity & Administration**
 - LAN / USB / Modem
- **On-instrument Data Storage and Logger**
- **Updated & Transferable calibrations**




**Grain
Flour
Food
Feed**

Automated Sample Presentation IM9140 Flour Analyzer

Benefits:

- **Extremely easy to use**
- **Consistent results – no operator influence**
- **Sample temperature measurement**
- **Reduces error & tampering**
- **No cups to pack**
 - Faster sample throughput
 - Improved accuracy & repeatability
 - No expensive replacement cups/cells



**Grain
Flour
Food
Feed**

DA 7200 Diode Array NIR

- **Breeders & Seed Companies**
 - High analysis throughput requirements
- **Use at Mills & Mix Plants**
 - Incoming ingredients
 - Optimize blends
 - Optimize additives
 - Monitor mixes
 - Moisture, Protein, Fat, Sugars
- **6 Second Analysis**
- **No sample preparation**
- **Non-contact measurement**
- **No cleaning**
- **Accurate, Flexible, Easy-to-use**



DA 7200 Application Examples

- **Whole Grain** – Moisture, Protein, Oil, Starch
- **Flour** – Moisture, Protein, Ash
- **Gluten** – Moisture, Protein
- **Midds** – Moisture, Protein
- **Bakery Mixes** – Moisture, Protein, Fat Sugars
(Pancake, Buns, Gravies, Brownies/Cakes)
- **Oils** – FFAs, Moisture, P, Linoleic Content
- **Butter** - Moisture, Fat, Salt
- **Icings / Fillings** – Moisture, Sugar, Fat
- **Cookies/Crackers** – Moisture, Fat, Fiber



**Grain
Flour
Food
Feed**

GIPSA USA Test, 9 Mills in 5 Countries

Quality Parameters	R ²	SECV	Range	Number of samples
Moisture, %	0.98	0.14	10.2-15.9	350
Flour protein, % db	0.96	0.26	8.3-15.6	350
Ash, % db	0.96	0.058	0.43-1.99	350
Ash, db < 0.8%	0.85	0.029	0.43-0.79	350
Color, Kent Jones Color grade	0.995	0.22	-3.5-4.8	140
Grain protein content, 14% mb	0.98	0.26	9.36-19.7	480
Flour protein content, 14% mb	0.99	0.18	8.24-17.84	350
SDS sedimentation vol., mL	0.68	2.55	24-45	270
Wet Gluten content %	0.92	0.32	16-34	270
Flour particle size GMD, 50% vol, µm	0.82	1.57	78.69-94.55	200
Insoluble glutenins, mg (IPP)	0.86	0.74	6.55-19.11	200
Total polymeric protein, mg (TPP)	0.94	0.63	45.88-63.81	200
Farinograph Absorption, %	0.80	1.5	55.5-73.1	250
Loaf volume, cm ³	0.90	34.2	685-1238	200
Loaf specific vol., cm ³ *g ⁻¹	0.87	0.24	4.62-8.24	200
Soluble glutenins, mg (mgSPP)	0.80	0.37	2.64-6.17	120
Soluble gliadins, mg (mgGLI)	0.91	0.54	6.81-15.7	120
Alveo. Length (L), mm	0.70	17.02	37-174	200
Alveo. Swelling index, mL	0.75	1.77	13.5-29.4	200
Starch Damage, UCD	0.94	3.94	-0.37-27.6	210
Starch Damage, AACC	0.93	0.35	1.06-6.66	135
Starch Damage, Ai	0.97	0.99	72.3-96.5	135

DA 7200 for Breakfast Cereals

Moisture content is critical for storage and crispiness.


Acceptable range is 2.7 – 3.7 %.

Rapid moisture meters are not accurate enough and oven measurement is too slow.



SimPLUS Software

- **Networking & Administration Software for the DA 7200**
- **Automatically update multiple instruments remotely**
- **Produce customized reports**
- **Integrate with LIMS systems and Formulation Software**
- **Results for other tests in the same report**
- **Action notification**



Certificate of Analysis

Sample ID: hylt Analyzed: 2007/10/16 5:26:18PM Address: Product: melle Customer: BASF Washington Supplier: Telephone: Line number:

Parameter	Unit	Result
1 Protein	%	8.3
2 Moisture	%	9.5
3 Fibre	%	2.6
4 Fat	%	0.6
5 Ca	%	0.02
6 P	%	0.22
7 lysine	%	0.22
8 Methionine	%	0.18
9 Cystine & Methionine	%	0.37
10 Threonine	%	0.29
11 Digestible Lysine Poultry	%	0.17

The Rapid Visco Analyzer (RVA)



An instrument designed for analysis of viscometric properties of materials under varying stirring (shear), heating, and cooling conditions


- **Rotational viscometer**
 - Stirs sample and measures resistance in centipoise (cP) units
- **Variable stirring speed**
 - Measure the sample as it will be processed
 - Add shear energy to see effect on the sample/matrix
- **Variable heating**
 - Measure the sample as it will be processed
 - Measure sample/matrix changes through temperature cycles
- **Stored profiles**
 - Allows samples to be analyzed under identical conditions




Milling & Baking: Flour Quality

- ICC 162, AACC 76-21
- 13 minutes pasting
- Bread, cake & noodle flours
- Measure sprouting, pasting potential, heat treatment, amylase addition
- Fingerprint flour premix


WHEAT FLOUR & MIX



Wheat Grain Soundness & Pasting



- Alpha-amylase activity
- Wheat, malting barley, rye



- 13 minutes pasting
- Different amylase activity


Rapid Visco Analyzers (RVA)

- **Wide variety of applications**
 - Starches, proteins, gums, hydrocolloids, in-process & finished products, process emulations
- **Flexible Set-up**
 - Shear, Temperature, Time
 - Automatically imitate complex plant processes
- **Miniature Pilot Plant**
 - Economical product development on a small scale
- **Ingredient Testing**
 - Monitor suppliers - See matrix effects before materials hit the plant floor In-Process Monitoring
 - Detect problems prior to finished products
- **Final Product Testing**
 - Satisfied customers
 - Product fingerprinting



doughLAB – Flour & Dough Rheology

- **New instrument for traditional rheological measurements**
 - Water Absorption, Mixing Tolerance, Stability, etc.
- **Easy Set-up & Cleanup**
- **Automated procedures with no glass**
- **doughLAB results are easy to read**
 - Curve analysis software increases accuracy & repeatability and removes operator dependence
- **doughMAP software provides calculations of blends to a specification automatically**
- **Greater flexibility than other dough mixers/testers**
- **New! – mimic a Farinograph**
 - Get Farinograph equivalent results from the doughLAB



What tests can doughLAB perform?

- **Mimic Standard Methods**
 - AACC 54-21 (300g flour, 14% moisture)
 - RACI 06-02
 - ICC 115/1
 - Constant speed (63rpm) & temperature (30°C)
- **Non-standard Methods**
 - Set up your own methods to fit your products and processes
- **Variable speed & temperature**
 - for R&D
 - imitate modern high energy mixing
 - imitate modern variable temperature mixing

**Grain
Flour
Food
Feed**


Micro-doughLAB

- Dough rheology on just 4 GRAMS



**Grain
Flour
Food
Feed**


Detecting bug damage with Gluten Index



Rosell et al.: J. Cereals Sci. 39 (2004) 187 - 193

**Grain
Flour
Food
Feed**

Gluten Index for Wheat and Flour



**Laboratory mill
LM 3100 or
LM 120**

**Glutomatic system:
Glutomatic 2200, Glutork 2020
& Centrifuge 2015**

**Grain
Flour
Food
Feed**

Visit us at

Perten
INSTRUMENTS

www.perten.com

 **NEWPORT**[®]
SCIENTIFIC

www.newport.com.au

Booth 44