



## **IAOM Meeting – Nairobi, 2018**

### New methods for raw material characterization



**Food Division | Food Application Lab | Alexandra Petz , Jessica Wiertz, Ulrike Ito**

**...where quality is measured.**

# New methods for raw material characterization

## Agenda

- Company profile
- Current food trends
- New methods




## Company profile

**...where quality is measured.**

# Brabender® GmbH & Co. KG

Part of the Brabender® Group



...where quality is measured.

**Brabender®**

**3** companies headquartered in Duisburg (Germany)

**100 %** family-owned

**6** subsidiaries abroad

**450** employees

**98** international sales agents

**80-85 %** export

Present in **129** countries.  
Steadily growing.



# Brabender® GmbH & Co. KG

## Portfolio

- Laboratory mills for sample preparation
- Instruments for grain reception tests (e. g. moisture testing)
- Instruments for flour quality tests
- Rheological equipment for dough quality testing
- Starch quality testing instruments
- Laboratory extruders



...where quality is measured.

**Brabender®**

# Brabender® GmbH & Co. KG

## Facilities

- Food Rheology Application Lab
- Food Extrusion Lab
- Bakery Lab



...where quality is measured.

**Brabender®**



## Current food trends

**...where quality is measured.**

# Current Food Trends

## Market trends & customer requirements

### Nutritional trends / technologies


- Gluten-free
- Low-carb
- Alternative protein sources
- Vegan food
- Ancient grains
- Whole meal products
- Superfood
- Algae
- Insects
- 3d-printed food

### Customer wishes

- Automatisation
- Quick methods
- Cost savings
- Independence from operator

### Conditions

- Legislation on fertilisation
- Hygienic requirements



## New methods

...where quality is measured.

...where quality is measured. **Brabender®**

# New methods

## Analysis of gluten-free raw materials

- Raw materials: corn flour, rice flour, flour made from pulses, etc.
- Initial situation: gluten-free raw materials cause doughs with non-elastic behaviour
- Solution: Farinograph-TS + Farino-Add S30  
Adaption of standard methods



...where quality is measured.

**Brabender®**

# New methods

## Analysis of gluten-free raw materials

### Wheat



www.shutterstock.com - 113273476

### Buckwheat






www.shutterstock.com - 715016584

### Soy



www.shutterstock.com - 1073226083



...where quality is measured.

**Brabender®**

# New methods

## Analysis of gluten-free raw materials

Gluten-free

Rice



www.shutterstock.com - 166476173

Brown rice






www.shutterstock.com - 375159193

Corn



www.shutterstock.com - 746658616



...where quality is measured.

**Brabender®**

# New methods

## Analysis of gluten-free raw materials

### Quinoa



www.shutterstock.com · 306455465

### Sorghum






www.shutterstock.com · 734380915

### Chickpea



www.shutterstock.com · 1031228197



...where quality is measured.

**Brabender®**

# New methods

## Analysis of gluten-free raw materials




**Red lentil**



**Coconut**



**Cricket**



...where quality is measured.


**Brabender®**

# New methods

## Analysis of gluten-free raw materials



### Corn flour + 2% HPMC



...where quality is measured.

**Brabender®**

# New methods

## Brabender GlutoPeak

- Raw materials: wheat flour, whole wheat flour, dried gluten
- Initial situation: requirements for quick methods at raw material reception
- Solution: GlutoPeak  
Adaption of standard methods



Method	
Amount of sample	9 g
Amount of liquid (dest. water)	9 g
Temperature	36 °C
Speed	2.750 rpm
Time	5 min

...where quality is measured.

**Brabender®**

# New methods

## Brabender GlutoPeak

### Principle:

- Rheologic measurement
- High-speed-mixing
- Torque measurement
- Gluten aggregation



### Target group:

- Grain breeders & traders
- Mills and bakeries
- Starch producers
- Research institutes

...where quality is measured.

**Brabender®**

# New methods

## Brabender GlutoPeak

### Which advantages does the GlutoPeak provide?

- Time saving due to quick methods
- Quick raw material control at grain or flour reception
- Limited dependence from operator due to easy handling
- Cost saving due to quick decisions regarding quality
- Easier silo management

*Brabender GlutoPeak  
successfully used by our reference customers*




...where quality is measured.

**Brabender®**

# New methods

## Analysis of dried vital gluten

- Project: **Development of a method to analyze dried vital gluten quality**
- Target: Being able to differentiate among qualities at raw material reception  
Having a measurement system independent from the operator
- Institution(s): Brabender, ttz Bremerhaven, CSM Deutschland GmbH, Blattmann Schweiz AG
- Result: RGC – Rapid Gluten Check\*



\*Gall, S. et. al. (2018)

...where quality is measured.


**Brabender®**

# New methods

## Analysis of dried vital gluten

### Method\*

Amount of sample (dried gluten)	2,1 g
Amount of liquid (dest. water)	4,4 g
Temperature	36 °C
Speed	500 rpm (1min) 0 rpm (2 min) 3.300 rpm (10 min)
Time	13 min



\*Gall, S. et. al. (2018)

...where quality is measured.

**Brabender®**

# New methods

## Analysis of dried vital gluten

### Target groups:

- Gluten producers
- Starch producers
- Gluten processors (bakery ingredients, baking mixtures, meat analouge producers)

### Application:

- Raw material quality check
- Final product check
- Process control during starch/gluten manufacturing
- Proof of heat damage

### Technical Note:<sup>\*</sup>

**Brabender**  
Qualität ist messbar.

**CSM**  
Bakery Solutions

**ttz Bremerhaven**

**Technical Note:**  
Development of a method to analyse vital wheat gluten with the GlutoPeak  
Sergej Gall<sup>1</sup>, Alexander Heinrich<sup>2</sup>, Julian Huen<sup>3</sup>, Jessica Wiertz<sup>1</sup>, Leslie Zeller<sup>1</sup>  
<sup>1</sup> Brabender GmbH & Co. KG, <sup>2</sup> CSM Deutschland GmbH, <sup>3</sup> ttz Bremerhaven

**GlutoPeak – new “Rapid Gluten Check” method**  
Vital wheat gluten is used as a powdered additive in bakery products to improve their quality, especially in order to achieve a consistent quality in the bake or to compensate for quality fluctuations in the flour.  
The properties of vital wheat gluten are significantly determined by its components, gliutenins and gliadins. While gliutenins create mechanical strength, the gliadins are responsible for the elastic properties (Schärf and Köller, 2016). The distribution of the gliutenins in relation to the gliadins, weight and the intergliutenin/gliadine bonds are an important criterion for determining dough and baking properties and the gluten quality (Ortolan and Steel, 2017). For this reason it has been very difficult and time-consuming until now to precisely determine the quality of the gluten. The following illustration shows the roles of both gliutenins and gliadins in the gluten network:




Fig. 1: Roles of gliutenins and gliadins in the vital wheat gluten network (Schärf and Köller, 2016)

**Challenges when analysing vital wheat gluten quality**  
The quality of powdered vital wheat gluten depends in turn on various factors. On the one hand, it is dependent on the quality of the cereal flour which it is produced. On the other hand, it is also dependent on various process parameters, for instance the extraction process of starch and gluten, the drying parameters, particularly time and temperature. However, the exact dependences have not yet been fully explained.  
Vital wheat gluten dimensions can often differ in terms of their quality. So far no standardised methods have been established to measure the quality of vital wheat gluten in dried form and to distinguish the product qualities of different batches. Baking trials give a good impression about quality, yet are very time-consuming to carry out. This poses a major challenge for the manufacturers and processors of vital wheat gluten as they have to be able to assess quality fluctuations and even them out.

\*Gall, S. et. al. (2018)

# New methods

## Analysis of whole wheat flours

- Project: **Development of a method to analyze whole wheat flour**
- Target: Having a quick method with saving time producing refined flour
- Institution(s): CWB, Wheat Marketing Centre, University of Minnesota, Washington State University, Anhui Agricultural University
- Result: Whole wheat method\*



\*Wang et. al. (2018)

...where quality is measured.

**Brabender®**

# New methods

## Analysis of whole wheat flours

### Method

Amount of sample (whole wheat flour) 8,0 g

Amount of liquid (CaCl<sub>2</sub> solution, concentration: 55,49 g/l) 10,0 g

Temperature 20 °C

Speed 3.000 rpm

### Effect of particle size:

- Particle size has got a significant effect on the analysis
- The method is more efficient when the particle size of the bran is smaller

\*Wang et. al. (2018)

...where quality is measured.

paper:<sup>\*</sup>



#### GlutoPeak method improvement for gluten aggregation measurement of whole wheat flour

Jungang Wang<sup>a,b</sup>, Gary G. Holt<sup>a,c</sup>, Ting Liu<sup>a,c</sup>, Naifu Wang<sup>a,d</sup>, Jayne Bock<sup>a</sup>

<sup>a</sup> Wheat Marketing Center, Inc., 1200 NW Market Parkway, Suite 230, Portland, OR 97209, USA

<sup>b</sup> Department of Biological Systems Engineering, Washington State University, Pullman, WA 99124, USA

<sup>c</sup> Department of Biological Sciences and Mathematics, Washington State University, Pullman, WA 99124, USA

<sup>d</sup> Department of Plant and Soil Sciences, Agricultural Experiment Station, Penn State, University Park, PA 16802, USA

<sup>e</sup> C.W. Brabender Instruments, Inc., 50 E. Wesley Street, South Hackensack, NJ 07606, USA

#### ARTICLE INFO

Keywords:  
Whole wheat flour  
Gluten  
Flour particle size  
Gluten aggregation  
Dough rheology

A shear-based device, the GlutoPeak, was developed to measure the aggregation behavior of gluten. In this study, the GlutoPeak testing method was improved and applied for the gluten aggregation evaluation of whole wheat flour samples over a wide range. The sample/water weight ratio, mixing time, and mixing temperature were adjusted to different levels to obtain the optimum results from the GlutoPeak. The Maximum Torque (MT) using two commercial WWF samples. The temperature had a greater effect on the improvement of the MT than the mixing time. The mixing time had a greater effect on the improvement of both PMT and MT. Another four commercial wheat samples were tested as WWFs of different particle size ranges. The gluten aggregation of the WWF was measured using the improved GlutoPeak parameter settings. The MT increased significantly with increasing temperature up to 30 °C. The PMT significantly increased with a decrease of WWF particle size. Compared to the MixLab and Farinograph tests, the GlutoPeak method took less time and provided significantly different ( $p < 0.05$ ) results for WWFs of different particle size ranges.

#### 1. Introduction

Whole wheat flour (WWF) is milled from the whole grain kernel that contains bran, germ, and endosperm, thus, it provides a highly nutritious and inexpensive source of protein, carbohydrates, dietary fiber, and various minerals. Considerable scientific research has been conducted on the properties of WWF (Kilmer & Holt, 2007; United States Department of Health & Human Services, 2007; Okaror & Liu, 2010). United States Department of Agriculture (USDA) recommends that Americans should consume at least half of their grains as whole grains (USDA, 2013), and most of the whole grain intake is whole wheat (USDA, 2013). Along with the increasing health awareness of consumers and recommendations of scientists, a large market demand for various WWF-based products is expected (Okaror & Liu, 2010).

One of the major challenges of WWF application is the influence of wheat bran, which results in changes in dough rheological properties and further affects the quality of food products (Cai, Cho, Hyun, Jeong, & Kim, 2011; Cho, Cai, Lee, Cho, Hyun, & Kim, 2014; Cho, Medina, Trippard, & Langton, 2015; Wong, Hsu, Kawan, & Lin, 2016). Many studies have been conducted to reduce the influence of bran in

products by adding vital wheat gluten, enzymes, and adjusting the particle size of bran (Kattus, Ishmukhamed-Murtila, Partana, Formell, & Antik, 2006; Pessina, Collier, & Harro, 2009). Compared to other methods, adjusting the particle size of bran is more natural and cost effective.

Several rheological instruments, such as the MixLab, Farinograph, and Mixograph, have been developed and applied to dough properties (Ghosh et al., 2010; Koksal, Kaderlik, Sarioglu, Ozel, & Dilek, 2009; Monthey & Scherzer, 2002; Vlachis & Danzik, 2013). However, it is difficult for these instruments to detect meaningful differences in dough rheological properties for WWFs of different particle sizes.

Previous studies reported no significant effects of bran particle size regarding water absorption in the Farinograph (Bauer, van Heege, & van der Steene, 1997) and the MixLab (van Heege, 1997) and dough development time in Mixograph (Cai et al., 2014) or MixLab (Lin, Hsu, Lee, Marques, & Dutier, 2010) of WWFs.

Recently, new instruments have been developed as a more efficient and accurate predictor of wheat bran functional quality (Caneil & Soethierman, 2012; Mart, Corchado, Frigola, Denisse, & Fagundes, 2014; Melchi, Denisse, Bouceta, Marcone, & Soethierman, 2011).

<sup>a</sup> Corresponding author. Wheat Marketing Center, Inc., 1200 NW Market Parkway, Suite 230, Portland, OR 97209, USA.

E-mail address: jwang@wmc.org (J.G. Wang).

https://doi.org/10.1016/j.lwt.2017.11.009

Received 2 October 2017; Received in revised form 30 November 2017; Accepted 30 November 2017

0020-0439/\$ - see front matter © 2017 Elsevier Ltd. All rights reserved.

**Brabender®**

# New methods

## Analysis of ancient grains

- Project: **Development of a method to analyze flour of ancient grains such as Einkorn, Emmer and Spelt**
- Target: Differentiation among standard wheat flours and ancient grain flours
- Institution(s): Leibniz-Institute for Food Systems Biology at the Technical University of Munich; Biotask AG
- Result:
  - Showing differences in aggregation behavior
  - Deriving conclusions on the baked product's quality




\*Geislitz, S. et. al. (2018)

# New methods


## Analysis of ancient grains

### Method\*

Amount of sample (flour)	8,5 g
Amount of liquid ( $\text{CaCl}_2$ solution, concentration: 55,49 g/l)	9,5 g
Temperature	34 °C
Speed	1.900 rpm
Time	7 min.



paper:\*\*



Gluten protein composition and aggregation properties as predictors for bread volume of common wheat, spelt, durum wheat, emmer and einkorn

Sabrina Geisslitz<sup>1</sup>, Herbert Wieser<sup>1</sup>, Katharina Anne Scherf<sup>2</sup>, Peter Koehler<sup>3,b,\*</sup>

<sup>1</sup>Leibniz Institute for Food Systems Biology at the Technical University of Munich, Lise-Meitner-Strasse 34, 85354, Freising, Germany

<sup>2</sup>Brabender AG, Schleißheimer Straße 14–16, 73758, Ellwangen an der Jagst, Germany

#### ARTICLE INFO

Keywords:  
Bread volume  
Emmer  
Einkorn  
Gliadins  
Wheat

#### ABSTRACT

The technological properties of the ancient wheat species spelt, emmer and einkorn are considered inferior compared to common wheat and durum wheat, whereas there are only few comparative studies between ancient and modern wheat species to help fill this gap. Protein content and aggregation properties of dough and bread properties were determined in a series of eight cultivars of common wheat, spelt, durum wheat, emmer and einkorn grown under standardized conditions at one location in the same year. Spearman correlation and principal component analysis (PCA) revealed that especially the content of gliadins, high molecular weight proteins and their aggregation properties were correlated with the bread volume of loaves made from wholemeal flours of the five wheat species using milled baking items. Based on their proximity to common wheat in the PCA, one cultivar each of spelt, emmer and einkorn was identified that had similar protein analytical, functional and rheological properties to common wheat. Furthermore, the relationship between aggregation behavior using the Gliadinsink test (GPT) enabled an estimation of dough properties and bread volume. Therefore, the fast and easy GPT may serve as an alternative to time consuming and labor intensive baking tests.

#### 1. Introduction

In the 21st century, the production area of the "ancient" (billed) wheat species einkorn (*Triticum monococcum* L., diploid), emmer (*T. dicoccum* L., tetraploid) and spelt (*T. spelta* L., hexaploid) is negligible compared to that of the modern (naked) wheat species common wheat (*T. aestivum* L., hexaploid) and durum wheat (*T. durum* L., tetraploid). Especially in the last decades ancient wheat species were replaced by modern wheat species due to higher grain yields (spelt 37%, emmer and einkorn 62% more yield compared to common wheat). Ancient species are still used with the exception of einkorn which has to be separated from the grain in the mill (Geislitz et al., 2013). Because some consumers associate the consumption of ancient wheats with health benefits, ancient wheat species have been bred for at least 10 years in the bread and brewing industry such as bread, pizza and beer have been developed (Geislitz et al., 2013). Studies on the contents of bioactive components (e.g. dietary fiber components, phenolic acids, folates) in ancient and modern wheats revealed only small differences in the contents of bioactive species. For example, spelt even though emmer and einkorn contained more of the celiac toxic gluten due to its yellow color (Elshay and Hey, 2013). More studies on a wider range of properties of ancient wheats grown under standardized conditions are currently needed to evaluate health benefits (Geislitz, 2013). Further advantages of ancient wheats include their disease tolerance, adaptation to different climatic conditions, low requirements of fertilizers and potential in genetic biodiversity (Lengin et al., 2013). In addition, traditional diploid wheats (einkorn) may contain lower amounts of immunogenic proteins and peptides compared to hexaploid species. For example, the celiac disease-active 33-mer peptide was not detected in emmer, durum wheat and einkorn wheats in contrast to about 10% of spelt, spelt wholemeal flour which had comparable contents of the 33-mer (Geislitz et al., 2013).

The baking quality of wheat flours is mostly determined by gluten quality and quantity. Gluten proteins are storage proteins and divided into two main groups: gliadins and glutelins. Gliadins (GA) are soluble in aqueous alcohol only after reduction of disulfide bonds.

Contents and composition of GLIA (α<sub>1</sub>, α<sub>1</sub>2,  $\alpha$ - and  $\gamma$ -GLIA) and GLUT (sub-gladiins, high- (HMW-GS) and low-molecular-weight-glutenins (LMW-GS)) are important parameters for bread volume. For example, spelt wholemeal flour had a higher protein content and a higher protein fractionation followed by reversed-phase high-performance liquid chromatography (RP-HPLC) (Wieser et al., 1998). One gliadin

\* Corresponding author.  
E-mail address: peter.koehler@lifsc.de (P. Koehler).

<sup>b</sup> Retired.

https://doi.org/10.1016/j.jcs.2018.04.012  
Received 7 June 2018; Received in revised form 23 July 2018; Accepted 21 August 2018  
Available online 27 August 2018  
0735-5210/© 2018 Elsevier Ltd. All rights reserved.

\*Marti, A. et.al (2015); \*\*Geisslitz, S. et. al. (2018)



...where quality is measured.

**Brabender®**

- Project: **Analyzing dough ingredients (e.g. sour dough, fat, salt or yeast) with the GlutoPeak**
- Target: To be able to make the effect of ingredients and their amount visible;  
To show the ingredients effect on the dough and the gluten aggregation
- Institution: University College Cork, Ireland
- Result\*: GlutoPeak can be used for product development
  - Definition of target curves for optimal receipts
  - Variation of single ingredients or their amount

\*Heitmann, M. et.al. (2017b)

### Effect of sour dough





- Lowering the pH-value
- Weakening the gluten network

\*Heitmann, M. et.al. (2017b)

...where quality is measured.

**Brabender®**

### Effect of salt and fat



- Salt strengthens the gluten
- Fat plasticized the dough

\*Heitmann, M. et.al. (2017b)

...where quality is measured.

**Brabender®**

# New methods

How food extruders can contribute to product development

## Extrusion & raw material analysis.... does that match??

Lab extruders for raw material control,  
product development and research

Torque measurement as Farinograph & Co.

Analysis of material properties under certain  
conditions (shear stress, temperature,  
pressure)



...where quality is measured.

**Brabender®**

# New methods

How food extruders can contribute to product development



Process	Application / product
Direct Expansion	Snacks, Cereals
Indirect Expansion	e.g. crab-crisps
Cold forming	Pellets, Drops
No Expansion	Pet Food, Feed, Pasta, Chewing Gum
Sweets	Licorice
Reactive Extrusion	Modified starch / flour




...where quality is measured.

**Brabender®**

# New methods

How food extruders can contribute to product development

## Texturized Vegetable Proteins (TVP)



...where quality is measured.

**Brabender®**


# New methods

## How food extruders can contribute to product development

### Texturized Vegetable Proteins (TVP)

Variation of product and process parameters such as:

- Moisture content
- Screw speed
- Screw configuration
- Temp. heating zones
- Temp. cooling die
- Diameter cooling die
- ...



...where quality is measured.

**Brabender®**

# New methods

How food extruders can contribute to product development



...where quality is measured.

**Brabender®**

# New methods

## How food extruders can contribute to product development

- Project: **Development of a procedure for expanded soy protein**
- Target: processing soy with the KE19 extruder
- Result: expanded product based on soy



3mm round die



...where quality is measured.

**Brabender®**

# ViscoQuick

## Brabender's newest viscometer



- Successor of the MVAG
- Universal viscosity measurement
- First implementation of new Brabender product design for table devices
- Precise, high reproducibility
- Compact, new architecture
- Flexible, different paddle geometries

...where quality is measured.


**Brabender®**

# ViscoQuick

## Measurement and Applications

### Applications:

- Starch
- Flour
- Backing goods
- Dough
- Sweets
- Milk
- Paper
- Cosmetics
- ...



...where quality is measured.

**Brabender®**

# Current Food Trends & Research Projects

## Future



...where quality is measured.

<https://3dprint.com/151348/barilla-pasta-3d-printer/>

**Brabender®**

# Current Food Trends & Research Projects

## Literature

- AACC (2018): Ancient Grains and Pseudocereals. <https://aaccpublications.aaccnet.org/doi/pdf/10.1094/CFW-61-5-0198> (Access at 1st September, 2018)
- Advanced Biotech (2018): The Ins and Outs of Cake and Cake Batter Flavours. [www.adv-bio.com/ins-outs-cake-cake-batter-flavors/](http://www.adv-bio.com/ins-outs-cake-cake-batter-flavors/) (Access at 1st September, 2018)
- Anusci, V. (2018): Food 3D Printing Starts from the Sweet Ending. <https://all3dp.com/food-3d-printing-starts-sweet-ending/> (Access at 1st September, 2018)
- Best, S. (2016): Is this the future of fine dining? Restaurant where all the food and even the table is 3D printed - but it will cost you £250 a head. <http://www.dailymail.co.uk/sciencetech/article-3685998/Is-future-fine-dining-Restaurant-food-table-3D-printed-cost-250-head.html> (Access at 25th August, 2018)
- Bernstein, M. (2018): The science of baking bread. <https://www.seriousseats.com/2014/10/breadmaking-101-the-science-of-baking-bread-and-how-to-do-it-right.html> (Access at 1st September, 2018)
- Bonduelle (2018): Legumino Pasta. <https://www.bonduelle.fr/nos-bons-produits/rayon-conserves-de-legumes/legumio-pasta--pois-chiches-et-mais-en-fusilli>; <https://www.bonduelle.fr/nos-bons-produits/rayon-conserves-de-legumes/legumio-pasta--lentilles-corail-et-carottes-en-penne> (Access at 25th August, 2018)
- Bundeszentrum für Ernährung (2018a): Proteintrend bei Lebensmitteln. <https://www.bzfe.de/inhalt/proteintrend-bei-lebensmitteln-32492.html> (Access at 25th August, 2018).
- Bundeszentrum für Ernährung (2018b): Pasta aus Kichererbsen. <https://www.bzfe.de/inhalt/pasta-aus-kichererbsen-co-32594.html>. (Access at 25th August, 2018).
- Bundeszentrum für Ernährung (2018b): Urgetreide. <https://www.bzfe.de/inhalt/urgetreide-28442.html> (Access at 25th August, 2018).
- Cullinley, C. (2014): Processing with ancient grains. <https://www.bakervandsnacks.com/Article/2014/01/23/Ancient-grains-fermentation-improvements-needed> (Access at 1st September, 2018)
- Eatbanza (2018): Collections. <https://www.eatbanza.com/collections/store?page=2> (Access at 25th August, 2018)
- EU Machinery (2018): Diosna Wendel. [http://eu-mach.ru/store/testomesy\\_promyshlennye/testomesilnye\\_mashiny\\_wendel/](http://eu-mach.ru/store/testomesy_promyshlennye/testomesilnye_mashiny_wendel/). (Access at 1st September, 2018)
- Halterman, T. (2014): Barilla Announces Their 3D Printed Pasta Contest Winners. <https://3dprint.com/32604/3d-printed-pasta-contest> (Access at 1st September, 2018)
- Low Carb Ernährung (2018): Die Low Carb Ernährung. <http://www.low-carb-ernaehrung.de/> (Access at 25th August, 2018)
- Menayang, A. (2016): **Tolerant Foods poised to be top contender in legume pasta category.** <https://www.foodnavigator-usa.com/Article/2016/11/03/Tolerant-Foods-poised-to-be-top-contender-in-legume-pasta-category> (Access at 1st September, 2018)
- Nestlé (2016): Forschen für den Fortschritt. <https://www.nestle.de/good2016/forschen-fortschritt> (Access at 25th August, 2018)
- Nieburg, O. (2013): Mixing ha major influence on gluten free bread. <https://www.bakervandsnacks.com/Article/2013/02/06/Gluten-free-bread-quality-with-mixing> (Access at 1st September, 2018)
- Scott, C. (2016): Barilla Continues to Develop Pasta 3D Printer, Envisions Gourmet Customization. <https://3dprint.com/151348/barilla-pasta-3d-printer/> (Access at 1st September, 2018)
- Thierry GmbH (2018): Viskosität. <https://www.thierry-qmbh.de/pruefgeraete/viscosity/> (Access at 1st September, 2018)
- Trenderous (2018): The 3d-printed food revolution is here and this is what it looks like. <http://www.trenderous.com/lifestyle/the-3d-printed-food-revolution-is-here-and-this-is-what-it-looks-like/>. (Access at 1st September, 2018)
- Waffel Meyer (2018): Industrie. <http://www.waffel-meyer.com/mobile/index.php?p=pindustrie> (Access at 1st September, 2018)
- Wholegrainsconcil (2018): Gluten free whole grains. <https://wholegrainscouncil.org/whole-grains-101/whats-whole-grain-refined-grain/gluten-free-whole-grains> (Access at 1st September, 2018)
- World Bakers (2017): Ingredient Report – From Ancient Grains and Pseudocereals to Modern Flour. <https://www.worldbakers.com/white-paper/ingredient-report-ancient-grains-pseudocereals-modern-flour/> (Access at 1st September, 2018)
- Yuan, L. (2015): WHEN WILL WE 3D PRINT OUR FOOD? A REPORT FROM THE WORLD'S FIRST 3D FOOD PRINTING CONFERENCE. <https://thisismold.com/process/cook/3d-food-printing-conference#.W4r9xUxuKuU> (Access at 1st September, 2018)

Thank you for your attention!

Visit us at our booth No. 68

– Your Brabender Food Sales & Application Team –



**Alexandra Petz**

**T: +49 (0) 203-7788-156**

**E: [Alexandra.Petz@Brabender.com](mailto:Alexandra.Petz@Brabender.com)**